

DC RECHARGEABLE FOR LABORATORY

M1



Controller	
Model	M1
Weight	185 g
Dimension	W 67 x D 132 x H 30(mm)
Input	AC 100 V ~AC 240 V
Output	DC 14.8 V , 1.2 A
Frequency	50 / 60 Hz
Charging time	About 2,5 hours
Using time	About 9 hours
Handpiece	
Model	M 110
Type	Brush
Speed	Up to 25,000 rpm
Weight	120 g
Dimension	Ø 24 x L 136(mm)

M1 E-type



Controller	
Model	M1
Weight	185 g
Dimension	W 67 x D 132 x H 30(mm)
Input	AC 100 V ~AC 240 V
Output	DC 14.8 V , 1.2 A
Frequency	50 / 60 Hz
Charging time	About 2,5 hours
Using time	About 9 hours
E-type motor	
Model	MAX-88ESP
Type	Brush
Speed	Up to 25,000 rpm
Weight	71 g
Dimension	Ø 22 x L 110(mm)
Angle, Handpiece	With all of ISO standard contra angle or handpiece