


# MICRO CAB<sup>+</sup>


Self Contained  
Dust Collection System



## NEW! EURO STYLE DUST CABINET.

**MicroCab+** offers the large internal work space and large hand guards of the larger MacroCab unit but with the smaller footprint of the MicroCab.

## THE MOST EFFECTIVE SANDBLAST AND GRINDING CONTAINMENT SYSTEM AVAILABLE.



**Wall Mountable!**  
Optional bracket allows wall mounting. Microcab is easily removed for cleaning.

REF # 91010 (120v)

REF # 91354 (220v)

- ◆ Increased efficiency, ultra quiet, fan/filter system.
- ◆ Easy to clean, pull-out filter system.
- ◆ Snap-in, disposable windows.
- ◆ Bright, fluorescent light.
- ◆ Built-in, low profile, carrying handle.
- ◆ Constructed of extruded aluminum and injection molded plastic
- ◆ Ultra compact: 8.7" (22 cm) wide x 10.5" (26.5 cm) deep x 7.9" (20 cm) high. Weight 7 lbs (3.2 kg).
- ◆ Input Voltage: 120V AC or 220V AC ~ 50/60Hz
- ◆ Power Consumption: 50 Watts

[www.danvillematerials.com](http://www.danvillematerials.com)

3420 Fostoria Way, Suite A-200 San Ramon, CA 94583 800.827.7940 925.973.0710 Fax 925.973.0764



Innovative Design. Reliable Performance

93295 Rev A